葡萄糖氧化酶法测定试剂盒 (Glucose Oxidase Method, GOD) E1010

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- 1、 Li D, Song J Z, Li H, et al. Storage lipid synthesis is necessary for autophagy induced by nitrogen starvation[J]. FEBS letters, 2015, 589(2): 269-276.
- 2、 Zhang J D, Han L, Yan S, et al. The non-metabolizable glucose analog D-glucal inhibits aflatoxin biosynthesis and promotes kojic acid production in *Aspergillus flavus*[J]. BMC microbiology, 2014, 14(1): 95.
- 3、 Yang Y, Li W, Liu Y, et al. Alpha-lipoic acid attenuates insulin resistance and improves glucose metabolism in high fat diet-fed mice[J]. Acta Pharmacologica Sinica, 2014.
- 4、 Wang L L, Qiao J, Liu H H, et al. Ratiometric Fluorescent Probe Based on Gold Nanoclusters and Alizarin Red-Boronic Acid for Monitoring Glucose in Brain Microdialysate[J]. Analytical chemistry, 2014, 86(19): 9758-9764.
- 5、 Huang L, Gong L, Jiang X, et al. Photoactivation of GLUT4 translocation promotes glucose uptake via PI3-K/Akt2 signaling in 3T3-L1 adipocytes[J]. Journal of Innovative Optical Health Sciences, 2014, 7(03).
- 6、 Yan P, Zhang L, Feng Y, et al. SHR3824, a novel selective inhibitor of renal sodium glucose cotransporter 2, exhibits antidiabetic efficacy in rodent models[J]. Acta Pharmacologica Sinica, 2014, 35(5): 613-624.
- 7、 Yang M, Dai J, Jia Y, et al. Overexpression of juxtaposed with another zinc finger gene 1 reduces proinflammatory cytokine release via inhibition of stress - activated protein kinases and nuclear factor - κB [J]. FEBS Journal, 2014, 281(14): 3193-3205.
- 8、 Zhang X, Yang R, Jia Y, et al. Hypermethylation of Sp1 binding site suppresses hypothalamic POMC in neonates and may contribute to metabolic disorders in adults: impact of maternal dietary CLAs[J]. Diabetes, 2014, 63(5): 1475-1487.
- 9、 Liu Y, He R. Fasting induces a high level of 3-nitrotyrosine in the brain of rats[J]. Neuroscience letters, 2010, 472(3): 204-209.
- 10、 Fu H, Pang S, Xue P, et al. High Levels of Expression of Fibroblast Growth Factor 21 in Transgenic Tobacco (*Nicotiana benthamiana*)[J]. Applied biochemistry and biotechnology, 2011, 165(2): 465-475.
- 11、 Zheng W, Feng X, Qiu L, et al. Identification of the antibiotic ionomycin as an unexpected peroxisome proliferator-activated receptor γ (PPAR γ) ligand with a unique binding mode and effective glucose-lowering activity in a mouse model of diabetes[J]. Diabetologia, 2013, 56(2): 401-411.
- 12、 Ma Y, Luo T, Xu W, et al. A new recombinant pituitary adenylate cyclase-activating peptide-derived peptide efficiently promotes glucose uptake and glucose-dependent insulin secretion[J]. Acta biochimica et biophysica Sinica, 2012: gms078.
- 13、 Zhang Z, Xue H L, Liu Y, et al. Yi-Qi-Zeng-Min-Tang, a Chinese medicine, ameliorates insulin resistance in type 2 diabetic rats[J]. World journal of gastroenterology: WJG, 2011, 17(8): 987.
- 14、 Yang Q, He Y, Wang W. The protective effect of Liu-Wei-Di-Huang-Fang in salt-sensitive hypertension rats[J]. Clinical and Experimental Hypertension, 2013, 36(6): 426-432.